# Plate & Screw System

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### 9) 수술법 안내(Surgical Technique Information):

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## ENGLISH

### 1) Device Description

The TDM Plates and Screw system consists of a family of flat and contoured plates and screws that make up the Mini and Mid Locking Plate and Screw System and the Small Locking Plate and Screw System. The Plates are constructed from Titanium alloy (Ti-6AL-4V) or pure Titanium (Ti) and come in a variety of configurations. The Plates are intended to be used with solid locking and non-locking screws and non-locking low profile Screws. The Screws are constructed from titanium alloy (Ti-6AL-4V) and are available as threaded locking screws, cortical or cancellous, Cannulated screw from 1.5mm to 7.5 in diameter and range from 6mm to 125 in length.

**Implant material specification:** Titanium alloy(Ti-6AL-4V) or Pure Titanium (Ti)

### 2) Indications

**Mini Mid Locking Plate and Screw System:**  
The Mini and Mid Locking Plate and Screw System is intended to be used in the hands, wrist, and small bones in the foot.

**Small Locking Plate and Screw System:**  
The small locking plate and screw system is indicated for the clavicle, scapula, olecranon, humerus, radius, ulna, tibia, and fibula, calcaneus and other small bones.

**Large Locking Plate and Screw System:**  
The large locking plate and screw system is indicated for the tibia, femur.

**Cannulated Screw System:**  
The Cannulated Screw System is intended to be used for fracture of small bones of the hand or foot (2.5mm) or small and large bones (4.0mm larger).

**Headless Compression Cannulated Screw System:**  
The Headless Compression Cannulated Screw System is intended to be used for a wide range of different indications in the hand, wrist and joint fusion (arthrodeses) in the foot (2.3mm & 3.5mm, 3.0mm & 4.0mm) and fixation of intra-articular fractures of the humerus, femur and tibia(3.5mm & 5.0mm).

**Limited Sliding Screw System:**  
The Limited Sliding Screw System is intended to be used for fracture fixation of the proximal femur, large bones and large bone fragments.

### 3) Contraindications:

- Do not use for surgeries other than those indicated.
- In case of material sensitivity, documented or suspected, appropriate tests should be performed for material suitability prior to implantation.
- Severe osteoporosis, compromised bone stock, insufficient or immature bone may not be suitable for use of this device.

- Any active or suspected latent infection, sepsis or marked local inflammation in or around the surgical area.
- Physical interference with other implants during implantation or use.
- Compromised vascularity, inadequate skin or neurovascular status.
- Patients who are unwilling or incapable of following post-operative care instructions.

#### 4) Adverse Effects:

The potential complications associated with the use of the TDM Plate and Screw system may include:

- Infections, both deep and superficial;
- Pain or discomfort;
- Foreign body reactions;
- Loosening, bending, cracking or fracture of implant components;
- Loss of anatomic position with nonunion or malunion may occur;
- Delayed correction in alignment; and
- Bone resorption or over-production.

#### 5) Warnings:

- The TDM implant is a single use device and must never be reused. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Only a trained physician with knowledge of orthopedic procedures should use this device. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
- The patient must be aware of temporary limitations of physical activity and danger of complications after procedure.
- In case of products that need bending during procedure, do not bend forcefully, try to keep the bending at a minimal level, and do not straighten the bent product since the product could be damaged.
- Post-operatively, the patient must be warned that failure to follow doctor's instructions may cause the following: delay in bone union caused by overweight, reoperation due to bone union failure or graft damage, etc.
- Metallic implants can loosen, fracture, corrode, migrate, cause pain, or stress shield bone even after a fracture is healed, particularly in young, active patients. The surgeon must make the final decision on implant removal if either of these occurs.
- Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended purpose, any of the following complications may occur:
  - (1) Corrosion and possible increased risk of infection;
  - (2) Migration of implant position resulting in injury;
  - (3) Risk of additional injury from postoperative trauma;
  - (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult;
  - (5) Pain, discomfort, or abnormal sensations due to the presence of the device; and
  - (6) Bone loss due to stress shielding.
- The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Any decision to avoid the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.

#### 6) Precautions:

- It is the responsibility of the physician to rely on his professional education, training and judgment to choose the most appropriate device and treatment option.
- This device is provided non-sterile and must be sterilized prior to use following the recommendations provided in this document.
- Inspect device for any noticeable physical deformations, cracks, pinholes, foreign substances, prior to use.

#### 7) Storage Conditions:

The TDM Plate and Screw system is provided non-sterile and should be stored in a clean, dry environment. The devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions over time.

#### 8) Sterilization and Cleaning:

The TDM Plate and Screw System is provided non-sterile. The implants are single use devices and need to be steam sterilized (autoclaved) in the tray provided by TDM using the conditions and parameters specified below.

Method	Sterilizer	Temperature	Exposure Time	Drying Time
Moist Heat Sterilization	Gravity	132°C (270°F)	15 minutes	30 minutes
	Pre-Vacuum	132°C (270°F)	4 minutes	30 minutes

The **accessory instruments** are reusable and must be cleaned after each use and re-sterilized (using the same parameters above) prior to each use. Recommended cleaning instructions for manual and automated cleaning are provided below:

**Note:** Clean instruments as soon as possible after use. If cleaning must be delayed, immerse instruments in a compatible detergent solution, spray with an instrument pre-soak solution, or cover instruments with a towel moistened with purified water to prevent drying and encrustation of surgical soil. It is recommended that instruments are reprocessed as soon as is reasonably practical following use.

**Note:** Use FDA-cleared sterilization wraps or another appropriate FDA-cleared accessory that has been validated to allow sterilant penetration and to subsequently maintain sterility.

#### Preparation for cleaning:

When properly performed, cleaning, disinfection and/or sterilization do not compromise the use and mechanical performance of these instruments. These instruments are used with or on patients who may harbor both recognized and unrecognized infections. To prevent the spread of infection, all reusable instruments must be thoroughly cleaned, disinfected, and sterilized after use on each patient.

1. No assembly/disassembly of these instruments is required unless stated on the labeling, directions for use, or literature assembly instructions (LAI) pertaining to cleaning, disinfection, and sterilization.
2. Devices that require disassembly should be disassembled prior to cleaning.
3. Remove dried-on soil from devices, especially in areas such as joints and crevices, prior to washing.

#### Manual Cleaning:

- 1) Prepare an enzymatic cleaning solution in accordance to the manufacturer's instructions.
- 2) Soak soiled devices for a minimum recommended time specified by the enzymatic cleaning solution manufacturer or 5 minutes, whichever is longer.
- 3) Prepare a pH neutral (pH 7-9) detergent cleaning solution in accordance to the manufacturer's instructions.
- 4) Use a soft non-metallic bristle brush to thoroughly scrub all traces of blood and debris from the device surfaces for one minute.

For All instruments except the instrument with lumen and the instrument with movable parts, move to Step 5.

#### For the instrument with Lumens

- Fully immerse the instrument in freshly prepared enzymatic or cleaning detergent solution to avoid aerosol generation.
- Use a long, tight fitting, soft lumen brush (plastic bristles, like nylon) to thoroughly scrub the lumen or cannula of each device. Minimum brush diameters are specified in Additional Technical Information for each lumen size.

- Push the brush through the entire length of the lumen using a twisting motion to remove debris.
- Clean the brush in the cleaning solution between each pass of the brush through the lumen by manually immersing the brush in the cleaning solution until no sign of debris can be seen.
- Pass the brush through the lumen as described above several times for 1 minute from both ends, cleaning the bristles before retracting until brush comes out clean.
- Use a syringe, or equivalent, filled to capacity with cleaning solution to flush the lumen of each device.

#### For the instruments with Movable Parts

- Fully immerse the instrument in freshly prepared enzymatic or cleaning detergent solution to avoid aerosol generation.
  - Thoroughly brush with a soft non-metallic bristle brush to remove all traces of blood and debris.
  - Pay close attention to threads, crevices, seams, and any hard to reach areas.
  - Actuate any moveable mechanisms, such as hinged joints, box locks, or spring-loaded features, to free trapped blood and debris.
  - If the components of the instrument can be retracted, retract or open the part while cleaning the exposed area.
  - Use a syringe, or equivalent, filled to capacity with enzymatic cleaning solution to flush the threads, crevices, seams, and any hard to reach areas.
- 5) Rinse the device with warm, (30°C-40°C), tap water for a minimum of one minute and until visual evidence of debris, soil, and cleaning solution are gone.
  - 6) Pay particular attention to thoroughly flush lumens, articulating areas, and flexible segments with warm, tap water.
  - 7) Ultrasonically clean the device components for 10 minutes in neutral pH detergent (pH 7-9), prepared in accordance with the manufacturer's instructions.

**NOTE:** Ultrasonic cleaning is only effective if the surface to be cleaned is immersed in the cleaning solution. Air pockets will decrease the efficacy of ultrasonic cleaning. Be sure to minimize air pocket or bubble formation by flushing lumens, cavities, crevices or springs with cleaning solution while the instrument is immersed in the ultrasonic cleaner tank.

- 8) Rinse the device components with warm, tap water for a minimum of one minute and until visual evidence of debris, soil, and cleaning solution are gone.
- 9) Perform a final rinse with Reverse Osmosis Deionized (RODI) or Purified (PUR) Water.
- 10) Dry the device components immediately after final rinse with a clean towel or clean compressed air until visibly dry

#### Automatic Washing:

- 1) Prepare the enzymatic detergent according to the manufacturer's recommendations:
  - Fully immerse the device in the prepared detergent and allow the device to soak for a minimum of five minutes. Actuate any moveable mechanisms, such as hinged joints, box locks, or spring loaded features, to free trapped blood and debris. If the components of the device can be retracted, retract or open that part whilst cleaning the area. For devices with flexible shafts, bend or flex the instrument under the cleaning solution while brushing the flexible areas.
  - Using a soft-bristled brush, scrub the device paying particular attention to crevices and other hard-to reach areas. Use a syringe to flush the holes and any difficult to reach areas.
  - Rinse the device with deionized water (DI) at ambient temperature to rinse the detergent. Actuate any movable parts while rinsing and flush the holes and difficult to reach areas using a syringe.
- 2) Transfer the device, disassembled and/or configured as specified by TDM, into mesh baskets:
  - Load instruments so that hinges are open and cannulations and holes can drain.
  - Place heavier instruments at the bottom of containers. Do not place heavy instruments on top of delicate instruments.
  - For instruments with concave surfaces, such as currettes, place instrument with the concave surface facing downward to facilitate draining.
  - Avoid contact between devices.

Select the instrument cycle and ensure the following set of cycle parameters are correctly programmed:

Phase	Recirculation time (min)	Water Temperature	Detergent type and Concentration (if applicable)
Pre-Wash	02:00	Cold Tap Water	N/A
Enzyme Wash	02:00	Hot Tap water	Neutral pH Enzymatic detergent
Wash	02:00	66.0°C (151°F) (set point)	Neutral pH Enzymatic detergent
Rinse 1	03:00	Hot tap Water	N/A
Thermal Rinse	01:00	82.2°C (180°F)	N/A
Purified Water	00 :10	Treated water	N/A
Dry Time	07:00	115°C (240°F)	N/A

- Start the instrument cycle.
- Check devices for visible soil. Repeat cleaning from 1 if soil is visible and inspect again.

#### Cleaning Inspection

- Visually inspect each device under normal lighting paying close attention to hard to reach areas.
- For difficult to view design features, apply 3% hydrogen peroxide. (Bubbling is indicative of the presence of blood).

**NOTE:** Rinse the instruments thoroughly with warm water following hydrogen peroxide testing.

- If visible soil remains, repeat cleaning procedure again.

#### Disinfection

- Devices must be terminally sterilized per the recommended parameters above prior to surgical use.

#### 9) Surgical Technique Information:

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. TDM provides detailed surgical techniques in electronic format on their website. You may also contact your TDM representative for a paper copy or onsite demonstration.

#### 10) MRI Compatibility:

The TDM Plate and Screw System has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in MR artifacts and patient injury.

#### 11) Customer Service and Complaint Handling:

- Contact your TDM representative or TDM Co. Ltd., directly for additional information and complaints regarding the quality, safety of this product.
- In case of product malfunction, patient's serious health issue, death, or side effect, must contact and notify the headquarters with as much information as possible (product name, product number, serial number, patient's condition, etc.) through fax or telephone.

- This product is a medical device, and, therefore, only persons in charge of handling medical devices can sell and use this product.
- This product is a single use product, and, therefore, it is prohibited to reuse.

Caution

Batch Code

Catalogue Number

Manufacturer

Rx Only

Federal law (United States) restricts this device for sale by or on the order of the physician

Do not reuse

Date of Manufacture

Consult Instructions for Use

Non Sterile