


An anatomical illustration of a human skeleton in a blue-tinted style. A white circle and line callout points from the shoulder joint area to the main title.

PROXIMAL HUMERUS

LOCKING PLATE & SCREW SYSTEM

SYSTEM OVERVIEW



PROXIMAL HUMERUS LOCKING PLATE & SCREW SYSTEM SYSTEM OVERVIEW

The advanced principle of TDM USA is to positively impact the healthcare system by partnering with medical professionals to define and deliver value-generating solutions. The Proximal Humerus Locking Plate and Screw System is a component of the extremity fixation product family that constitutes the **ADVANCED PRINCIPLE IMPLANT SYSTEMS**.



THE PROXIMAL HUMERUS LOCKING PLATE AND SCREW SYSTEM offers multiple fixation options in a variety of configurations. It is intended for use in the proximal humerus and humeral shaft. The plates are constructed of titanium with a low profile proximal design and are used in conjunction with solid locking and non-locking screws. The screws are titanium alloy, available in diameters of 3.5mm and 4.0mm, and range in length from 20mm-60mm.

PROXIMAL HUMERUS LOCKING PLATE & SCREW SYSTEM

IMPLANT ORDERING INFORMATION

PROXIMAL HUMERUS LOCKING PLATE (PHLP)

PHLP 3H	83mm	194-10003
PHLP 4H	95mm	194-10004
PHLP 5H	107mm	194-10005
PHLP 6H	119mm	194-10006
PHLP 8H	143mm	194-10008
PHLP 10H	167mm	194-10010
PHLP 11H	179mm	194-10011



PROXIMAL HUMERUS LOCKING PLATE (PHLP), SHAFT

PHLP, Shaft 6H	78mm	194-10306
PHLP, Shaft 7H	90mm	194-10307
PHLP, Shaft 8H	102mm	194-10308
PHLP, Shaft 9H	114mm	194-10309
PHLP, Shaft 10H	126mm	194-10310
PHLP, Shaft 12H	150mm	194-10312
PHLP, Shaft 14H	174mm	194-10314



3.5mm CORTICAL SCREW

3.5mm Cortical Screw	20mm	235-21020
3.5mm Cortical Screw	22mm	235-21022
3.5mm Cortical Screw	24mm	235-21024
3.5mm Cortical Screw	26mm	235-21026
3.5mm Cortical Screw	28mm	235-21028
3.5mm Cortical Screw	30mm	235-21030
3.5mm Cortical Screw	32mm	235-21032
3.5mm Cortical Screw	34mm	235-21034
3.5mm Cortical Screw	36mm	235-21036
3.5mm Cortical Screw	38mm	235-21038
3.5mm Cortical Screw	40mm	235-21040



4.0mm LOCKING SCREW

4.0mm Locking Screw	20mm	240-29020
4.0mm Locking Screw	22mm	240-29022
4.0mm Locking Screw	24mm	240-29024
4.0mm Locking Screw	26mm	240-29026
4.0mm Locking Screw	28mm	240-29028
4.0mm Locking Screw	30mm	240-29030
4.0mm Locking Screw	32mm	240-29032
4.0mm Locking Screw	34mm	240-29034
4.0mm Locking Screw	36mm	240-29036
4.0mm Locking Screw	38mm	240-29038
4.0mm Locking Screw	40mm	240-29040
4.0mm Locking Screw	42mm	240-29042
4.0mm Locking Screw	44mm	240-29044
4.0mm Locking Screw	46mm	240-29046
4.0mm Locking Screw	48mm	240-29048
4.0mm Locking Screw	50mm	240-29050
4.0mm Locking Screw	52mm	240-29052
4.0mm Locking Screw	54mm	240-29054
4.0mm Locking Screw	56mm	240-29056
4.0mm Locking Screw	58mm	240-29058
4.0mm Locking Screw	60mm	240-29060



PROXIMAL HUMERUS LOCKING PLATE & SCREW SYSTEM

INSTRUMENT OVERVIEW

Drill Bit

2.5 X 130mm (White) 901-01125
3.5 X 160mm (Green) 901-22135



Drill Bit, AO QC

2.5 X 130mm (White) 901-31025
3.5 X 160mm (Green) 901-31135



Locking Drill Sleeve

3.5mm (Green) 901-12235



Locking Drill Sleeve, Single Handle

3.5mm (Green) 901-02235



Drill Sleeve, Single Handle

2.5mm (White) 901-13325
3.5mm (Green) 901-13335



DCP Drill Guide

2.5mm 901-03035



Drill Guide Block

901-10463



Depth Gauge

2.7/3.5/4.0mm 901-04027



Screwdriver, Self-Retaining

2.5mm Hex 901-35025



Screwdriver Shaft, Self-Retaining, AO QC

2.5mm Hex 901-15425



Torque Limiting Handle, AO QC

2.0Nm 901-17020



T-Handle, AO QC

910-06017



Forceps

901-08001

