



# Multi MIS Retractor System

# MULTA MIS RETRACTOR SYSTEM

The One Surgical **Multa MIS Retractor System** is a pedicle-based lumbar access system designed to support MIS TLIF procedures. This retractor system provides full adjustment control, making it highly versatile for MIS fusion procedures.

MAIN FRAME ADJUSTS TO  
PATIENT DEPTH



55 - 95mm  
LATERAL BLADE  
LENGTHS

NARROW AND  
WIDE LATERAL  
BLADES



TRANSLATION AND DISTAL  
TOE ON ALL ARMS

USED IN CONJUNCTION WITH THE  
MULTA MIS PEDICLE SCREW SYSTEM