



Magna L IBF System

SPINE

MAGNA L IBF SYSTEM

The One Surgical **Magna L Interbody Fusion System** features a state-of-the-art, engineered, self-distracting bulleted nose for optimal insertion. The anatomically-shaped endplate interface maximizes stability in PLIF and TLIF procedures. Each implant is precision machined from either specially-formulated PEEK polymer or Titanium, which incorporates surface texturing intended to promote bone on-growth.

TLIF


5° LORDOTIC INTERFACE

LARGE AXIAL
GRAFT WINDOW

TAPERED, SELF-DISTRACTING
BULLET NOSE

PLIF

BICONVEX,
7° LORDOTIC INTERFACE



6 - 16mm HEIGHTS*
22, 28, AND 32mm LENGTHS



6 - 16mm HEIGHTS*
22, 26, AND 30mm LENGTHS

*1mm increments