



# Multi Pedicle Screw System

# MULTA PEDICLE SCREW SYSTEM

The One Surgical **Multa Pedicle Screw System** is designed to fully lock at 65 in-lbs (7.3 Nm), nearly half the torque required by most competitive designs. The low torque design also fosters high tactile feedback and reduces anatomical loads during final construct tightening. The unique thread form incorporates a progressive, multi-diameter design intended to maintain purchase with changes in screw position.



UP TO 60° POLYAXIAL ANGLULATION



4.5 - 8.5mm SCREW DIAMETERS  
25 - 100mm LENGTHS

