

Clavicle Locking Plate & Screw System

CLAVICLE LOCKING PLATE & SCREW SYSTEM

The One Surgical **Clavicle Locking Plate & Screw System** offers multiple fixation options in a variety of configurations. The low profile plates are constructed of titanium alloy and are used with locking and non-locking screws, ranging in length from 10-30mm. The screws are available in 3.5mm diameter for shaft placement and 2.7mm diameter for distal placement.


CLAVICLE LOCKING PLATE OPTIONS

 <p>CLAVICLE LOCKING PLATE, MID SHAFT 4 to 9 Hole Left and Right</p>	 <p>CLAVICLE LOCKING PLATE, ANTERIOR 1 to 8 Hole</p>
 <p>CLAVICLE LOCKING PLATE, DISTAL 4 to 7 Hole Left and Right</p>	 <p>CLAVICLE LOCKING PLATE, DISTAL HOOK 2 to 4 Hole Left and Right</p>

2.7mm LOCKING SCREW
10 - 30mm (2mm increments)



2.7mm CORTICAL SCREW
Low Profile Head
10 - 30mm (2mm increments)



3.5mm LOCKING SCREW
10 - 30mm (2mm increments)



3.5mm CORTICAL SCREW
10 - 30mm (2mm increments)

