




CANNULATED/HCC SCREW SYSTEMS

SYSTEM OVERVIEW



CANNULATED/HCC SCREW SYSTEMS SYSTEM OVERVIEW

The advanced principle of TDM USA is to positively impact the healthcare system by partnering with medical professionals to define and deliver value-generating solutions. The Cannulated and Headless Compression Cannulated Screw Systems are components of the extremity fixation product family that constitutes the **ADVANCED PRINCIPLE IMPLANT SYSTEMS**.



THE CANNULATED SCREW SYSTEMS offer multiple sizes of standard and headless screws for a variety of applications. The titanium alloy screws address indications in the hand and foot, in small fragment applications, and in hip fractures.

CANNULATED/HCC SCREW SYSTEMS

IMPLANT ORDERING INFORMATION

2.5mm CANNULATED SCREWS

Short Threads (~1/3)

10 - 30mm (2mm increments)

725-200xx*



5mm Washer

291-20005



4.0mm CANNULATED SCREWS

Short Threads (~1/3)

20 - 40mm (2mm increments)

740-200xx*

45 - 70mm (5mm increments)



7mm Washer

291-20007



4.5mm CANNULATED SCREWS

Short Threads (~1/3)

20 - 50mm (2mm increments)

745-200xx*

54 - 78mm (4mm increments)



10mm Washer

(5.0mm Inner Diameter)

291-20110



5.0mm CANNULATED SCREWS

Short Threads (~1/3)

30 - 80mm (5mm increments)

750-200xx*



10mm Washer

(5.4mm Inner Diameter)

291-20010



6.5mm CANNULATED SCREWS

Short Threads (~1/3)

45 - 120mm (5mm increments)

765-20xxx*



12mm Washer

291-20012



2.3mm HEADLESS COMPRESSION CANNULATED SCREWS

10 - 30mm (2mm increments)

723-200xx*



3.0mm HEADLESS COMPRESSION CANNULATED SCREWS

10 - 30mm (2mm increments)

730-200xx*



3.5mm HEADLESS COMPRESSION CANNULATED SCREWS

16 - 40mm (2mm increments)




















735-200xx*



*xx = Screw length





















CANNULATED/HCC SCREW SYSTEMS

INSTRUMENT OVERVIEW

2.5mm Cannulated Screws				
Guide Wire	1.0 x 130mm	904-30001		
Guide Wire Sleeve	1.0mm	904-30202		
Cannulated Drill Bit	1.7 x 90mm	904-30108		
Cannulated Depth Gauge	1.0mm	904-30003		
Cannulated Screwdriver	2.0mm Hex	904-31207		
Small Cannulated Chuck, AO QC		904-40012		
Cannulated Screwdriver Shaft, AO QC	2.0mm Hex	904-30107		
Screwdriver Shaft, AO QC	2.0mm Hex	904-30009		
Forceps		901-08001		
4.0mm Cannulated Screws				
Guide Wire	1.5 x 147mm	904-40001		
Guide Wire Sleeve	1.5mm	904-40202		
Cannulated Drill Bit	2.6 x 103mm	904-40203		
Cannulated Tap with T-Handle	4.0mm	904-40204		
Cannulated Countersink with Depth Gauge	4.0mm	904-45005		
Cannulated Screwdriver	2.5mm Hex	904-41107		
Small Cannulated Chuck, AO QC		904-40012		
Cannulated Screwdriver Shaft, AO QC	2.5mm Hex	904-41007		
Screwdriver Shaft, AO QC	2.5mm Hex	901-15425		
Forceps		901-08001		

















CANNULATED/HCC SCREW SYSTEMS

INSTRUMENT OVERVIEW

4.5mm Cannulated Screws				
Guide Wire	1.5 x 147mm	904-40001		
Guide Wire Sleeve	1.5mm	904-40202		
Cannulated Drill Bit	3.2 x 138mm	904-45032		
Cannulated Tap with T-Handle	4.5mm	904-45004		
Cannulated Countersink with Depth Gauge	4.5mm	904-45005		
Cannulated Screwdriver	3.5mm Hex	904-45007		
Small Cannulated Chuck, AO QC		904-40012		
Cannulated Screwdriver Shaft, AO QC	3.5mm Hex	904-45008		
Screwdriver Shaft, AO QC	3.5mm Hex	904-45010		
Forceps		901-08001		
5.0mm Cannulated Screws				
Guide Wire	2.0 x 163mm	904-50001		
Guide Wire Sleeve	2.0mm	904-50202		
Parallel Drill Guide	2.0mm	904-51015		
Cannulated Drill Bit	3.6 x 138mm	904-50103		
Cannulated Tap with T-Handle	5.0mm	904-50204		
Cannulated Countersink with Depth Gauge	5.0mm	904-50205		
Cannulated Screwdriver	3.5mm Hex	904-50107		
Cannulated Chuck, Large QC		904-50012		
Cannulated Screwdriver Shaft, Large QC	3.5mm Hex	904-50108		
Screwdriver Shaft, Large QC	3.5mm Hex	904-50110		
Forceps		901-08001		





















CANNULATED/HCC SCREW SYSTEMS

INSTRUMENT OVERVIEW

6.5mm Cannulated Screws				
Guide Wire	3.0 x 255mm	904-65001		
Multiple Guide Wire Sleeve	3.0mm	904-66002		
Parallel Drill Guide	3.0mm	904-66012		
Cannulated Drill Bit	5.0 x 230mm	904-65003		
Cannulated Tap with T-Handle	6.5mm	904-65104		
Cannulated Countersink with Depth Gauge	6.5mm	904-65205		
Cannulated Screwdriver	5.0mm Hex	904-65107		
Cannulated Chuck, Large QC		904-50012		
Cannulated Screwdriver Shaft, Large QC	5.0mm Hex	904-65108		
Screwdriver Shaft, Large QC	5.0mm Hex	904-65010		
Forceps		901-08001		
2.3mm Headless Compression Cannulated Screws				
Guide Wire	1.0 x 130mm	904-30001		
Guide Wire Sleeve	1.0mm	904-30202		
Cannulated Drill Bit	1.7 x 90mm	904-30108		
Cannulated Drill Bit	2.6 x 103mm	904-40203		
Cannulated Depth Gauge	1.0mm	904-30003		
Cannulated Screwdriver	2.0mm Hex	904-31207		
Small Cannulated Chuck, AO QC		904-40012		
Cannulated Screwdriver Shaft, AO QC	2.0mm Hex	904-30107		
Screwdriver Shaft, AO QC	2.0mm Hex	904-30009		
Forceps		901-08001		

CANNULATED/HCC SCREW SYSTEMS

INSTRUMENT OVERVIEW

3.0mm Headless Compression Cannulated Screws				
Guide Wire	1.0 x 130mm	904-30001		
Guide Wire Sleeve	1.0mm	904-30202		
Cannulated Drill Bit	1.7 x 90mm	904-30108		
Cannulated Drill Bit	3.2 x 138mm	904-45032		
Cannulated Depth Gauge	1.0mm	904-30003		
Cannulated Screwdriver	2.0mm Hex	904-31207		
Small Cannulated Chuck, AO QC		904-40012		
Cannulated Screwdriver Shaft, AO QC	2.0mm Hex	904-30107		
Screwdriver Shaft, AO QC	2.0mm Hex	904-30009		
Forceps		901-08001		
3.5mm Headless Compression Cannulated Screws				
Guide Wire	1.5 x 147mm	904-40001		
Guide Wire Sleeve	1.5mm	904-40202		
Cannulated Drill Bit	2.6 x 103mm	904-40203		
Cannulated Drill Bit	3.6 x 138mm	904-50103		
Cannulated Depth Gauge	1.5mm	904-40106		
Cannulated Screwdriver	2.5mm Hex	904-41107		
Small Cannulated Chuck, AO QC		904-40012		
Cannulated Screwdriver Shaft, AO QC	2.5mm Hex	904-41007		
Screwdriver Shaft, AO QC	2.5mm Hex	904-15425		
Forceps		901-08001		

EVERY PATIENT IS DIFFERENT EVERY SURGICAL FACILITY IS DIFFERENT

TDM USA offers multiple plate and screw fixation systems for orthopedic reconstruction and traumatic injuries of the upper and lower extremities.

And we partner with medical professionals to define and deliver value-generating solutions for today's complex healthcare challenges.

PRODUCT SOLUTIONS

