



Mini Locking Plate & Screw System

MINI LOCKING PLATE & SCREW SYSTEM

The One Surgical **Mini Locking Plate & Screw System** offers multiple fixation options in a variety of configurations. Plates in this system are intended for use in the small bones of the hand, wrist, and foot. The low profile plates are constructed of titanium and are used with solid locking and non-locking titanium alloy screws.



MINI LOCKING X-PLATE (5 Hole)
DOUBLE X-PLATE (8 Hole)
TRIPLE X-PLATE (11 Hole)



MINI LOCKING ZIG ZAG PLATE
5 and 7 Hole



MINI LOCKING DCP PLATE
4 to 8 Hole



MINI LOCKING RECONSTRUCTION PLATE
4 to 8 Hole



MINI LOCKING L-PLATE
3 to 8 Hole
Standard and Oblique
Left and Right



MINI LOCKING 3 T-PLATE
3 to 8 Hole



MINI LOCKING Y PLATE
3 to 8 Hole



MINI LOCKING CONDYLAR PLATE, PIN
3, 4, and 6 Hole
Left and Right



MINI LOCKING Y PLATE, WEDGE
2 and 3 Hole
3 - 6mm

1.5mm LOCKING SCREW
6 - 18mm (1mm increments),
20mm



2.0mm LOCKING SCREW
6 - 18mm (1mm increments),
20mm



1.5mm CORTICAL SCREW
6 - 18mm (1mm increments),
20mm



2.0mm CORTICAL SCREW
6 - 18mm (1mm increments),
20mm

