





Geminus Pedicle Screw System

GEMINUS PEDICLE SCREW SYSTEM


The One Surgical **Geminus Pedicle Screw System** is a low profile system featuring robust locking mechanisms and versatile applications. Intuitive tray layout simplifies surgical flow and speeds implantation. Over-molded instrument handles provide a comfortable surgeon experience.



CANNULATED VERSION ALSO AVAILABLE



Cortical Cancellous Screw



Dual Thread Screw

HYDROXYAPATITE COATING PROMOTES OSSEOINTEGRATION



Cortical Cancellous Screw - HA Coating



Cortical Cancellous Screw - HA Coating

HA coating adds .18mm to major diameter and .25mm to minor diameter